KBIS

- Upon selecting this program, we (the experimenters) should be able to specify:
 - 1. How long each screen will be presented for, in milliseconds; we will need to have control over each screen individually
 - 2. Whether the length of the line will be held constant from screen 1 to screen 3 (lower priority)
 - 3. The number of trials that will run successively
- For each trial, the program should record:
 - 1. targetNum: The location of the hashmark relative to the left end of the line presented in screen 1 (pixels is fine)
 - 2. lineLength: The length of the line presented in screen 3 (pixels is fine; recording this is actually only necessary when the length of the line changes from screen 1 to screen 3)
 - 3. touch: The location of the participant's touch relative to the left end of the line presented in screen 3 (pixels is fine)
 - 4. trialNum: The location of the trial within the sequence of trials run successively (1, 2, 3, etc.)
 - 5. stimScreen, delay1, respScreen, dispScreen, delay2: How long screen 1, screen 2, screen 3, screen 5, and screen 6 were presented for, respectively (ms)
 - 6. lengthChange: A binary variable indicating whether the length of the line was held constant from screen 1 to screen 3
 - 7. blockLength: The number of trials run successively (should be the same as the last "trialNum")
 - 8. reactionTime: The time between when screen 3 is displayed, and when the participant touches the screen (ms)
 - 9. retry: A binary variable indicating whether the participant retried a given trial
- If the length of the line presented in screen 1 is 1000 units:
 - 1. The possible locations of the hashmarks in screen 1 are: [25, 50, 75, 100, 125, 150, 175, 200, 225, 250, 275, 300, 325, 350, 375, 400, 425, 450, 475, 500, 525, 550, 575, 600, 625, 650, 675, 700, 725, 750, 775, 800, 825, 850, 875, 900, 950, 975]
 - 2. The possible lengths of the line in screen 3b are: [400, 450, 500, 550, 600, 900, 950, 1000, 1050, 1100, 1400, 1450, 1500, 1550, 1600]
- Other notes:
 - 1. The length of the line on screen 1 should be about half of the tablet's screen width
 - 2. The mark on the line is recorded when the participant first touches the screen, not when they release their finger
 - 3. If the participant does not respond quickly enough on a given trial, the experiment should advance to the next trial in the block, and record that the participant did not respond
 - 4. If the participant retries a trial and misses a second time, move to the next trial and record as missing data

